


BioVendor

Research
and Diagnostic Products



HUMAN VASPIN ELISA

Product Data Sheet

Cat. No.: RD191097200R

For Research Use Only

CONTENTS

| | | |
|-----|--|----|
| 1. | INTENDED USE | 3 |
| 2. | STORAGE, EXPIRATION | 3 |
| 3. | INTRODUCTION | 4 |
| 4. | TEST PRINCIPLE | 5 |
| 5. | PRECAUTIONS | 5 |
| 6. | TECHNICAL HINTS | 6 |
| 7. | REAGENT SUPPLIED | 6 |
| 8. | MATERIAL REQUIRED BUT NOT SUPPLIED | 7 |
| 9. | PREPARATION OF REAGENTS | 7 |
| 10. | PREPARATION OF SAMPLES | 9 |
| 11. | ASSAY PROCEDURE | 10 |
| 12. | CALCULATIONS | 11 |
| 13. | PERFORMANCE CHARACTERISTICS | 13 |
| 14. | DEFINITION OF THE STANDARD | 17 |
| 15. | PRELIMINARY POPULATION AND CLINICAL DATA | 17 |
| 16. | METHOD COMPARISON | 18 |
| 17. | TROUBLESHOOTING AND FAQs | 19 |
| 18. | REFERENCES | 20 |
| 19. | EXPLANATION OF SYMBOLS | 21 |

**»» This kit is manufactured by:
BioVendor – Laboratorní medicína a.s.**

»» Use only the current version of Product Data Sheet enclosed with the kit!

1. INTENDED USE

The RD191097200R Human Vaspin ELISA is a sandwich enzyme immunoassay for the quantitative measurement of human vaspin (visceral adipose tissue-derived serpin).

»» Features

- **It is intended for research use only**
- The total assay time is less than 3.5 hours
- The kit measures vaspin in serum and heparin plasma
- Assay format is 96 wells
- Standard is recombinant protein based
- Quality Controls are human serum based
- Components of the kit are provided ready to use, concentrated or lyophilized

2. STORAGE, EXPIRATION

Store the complete kit at 2-8°C. Under these conditions, the kit is stable until the expiration date (see label on the box).

For stability of opened reagents see Chapter 9.

3. INTRODUCTION

Visceral adipose tissue – derived serpin A12 (vaspin), also named OL-64, an adipocytokine, is structurally a member of the serine protease family. Serpins are the most diverse family of protease inhibitors. Their typical structural feature is the core domain composed from 3 beta-sheets and 9 alpha-helices. The inhibitory activity of vaspin has not been described up to now, but its reactive site loop is typical for this proteinase family. Human vaspin protein is composed of 395 amino acids and has a molecular weight of approximately 45.2 kDa and predicted pI 9.26.

The cDNA was first isolated from white adipose tissue of Otsuka Long-Evans Tokushima Fatty (OLETF) rats. Vaspin mRNA expression is specific for visceral adipose tissues and it is also found circulating in the serum. The level of serum vaspin increased with age up to the peak of obesity, body weight and insulin resistance in OLETF rats and decreases with worsening of diabetes. Vaspin expression is missing in the diabetes-resistant lean rats, LETO, in comparison to OLETF rats, animal model of metabolic syndrome. Expression was also absent in the subdermal, brown fatty tissue and other non-adipose tissues in OLETF rats. These findings lead to the conclusion that the target tissue for insulin sensitising effect of vaspin is white adipose tissue.

In humans, elevated serum concentration of vaspin is associated with obesity and impaired insulin sensitivity. In patients with type 2 diabetes the correlation between increased vaspin levels and BMI and decreased insulin sensitivity has not been observed. Vaspin expression decreased when diabetes worsened and its levels normalised when insulin or pioglitazone was administered.

Gender differences in vaspin serum levels have been found in separate studies by two different authors. The low levels of vaspin in serum seem to be typical for lean subjects and athletes with long-term physical training. On the other hand, serum vaspin concentration increased in overweight people after they lost weight because of increased exercise. This paradox has been explained by the fact that serum vaspin level is differentially regulated in the non-active resting state and after exercise.

Areas of investigation:

Metabolic syndrome

Obesity

Insulin resistance

4. TEST PRINCIPLE

In the Biovendor Human Vaspin ELISA, standards, quality controls and samples are incubated in microtitration wells pre-coated with polyclonal anti-human vaspin antibody. After a 60 minute incubation followed by washing, biotin labelled polyclonal anti-human vaspin antibody is added and incubated with the captured vaspin for 60 minutes. After another washing, streptavidin-HRP conjugate is added. After 30 minutes incubation and the last washing step, the remaining conjugate is allowed to react with the substrate solution (TMB). The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured. The absorbance is proportional to the concentration of vaspin. A standard curve is constructed by plotting absorbance values against vaspin concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

5. PRECAUTIONS

- **For professional use only**
- Wear gloves and laboratory coats when handling immunodiagnostic materials
- Do not drink, eat or smoke in the areas where immunodiagnostic materials are being handled
- This kit contains components of human origin. These materials were found non-reactive for HBsAg, HCV antibody and for HIV 1/2 antigen and antibody. These materials should be handled as potentially infectious, as no test can guarantee the complete absence of infectious agents
- This kit contains components of animal origin. These materials should be handled as potentially infectious
- Avoid contact with the acidic Stop Solution and Substrate Solution, which contains hydrogen peroxide and tetramethylbenzidine (TMB). Wear gloves and eye and clothing protection when handling these reagents. Stop and/or Substrate Solutions may cause skin/eyes irritation. In case of contact with the Stop Solution and the Substrate Solution wash skin/eyes thoroughly with water and seek medical attention, when necessary
- The materials must not be pipetted by mouth

6. TECHNICAL HINTS

- Reagents with different lot numbers should not be mixed
- Use thoroughly clean glassware
- Use deionized (distilled) water, stored in clean containers
- Avoid any contamination among samples and reagents. For this purpose, disposable tips should be used for each sample and reagent
- Substrate Solution should remain colourless until added to the plate. Keep Substrate Solution protected from light
- Stop Solution should remain colourless until added to the plate. The colour developed in the wells will turn from blue to yellow immediately after the addition of the Stop Solution. Wells that are green in colour indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution
- Dispose of consumable materials and unused contents in accordance with applicable national regulatory requirements

7. REAGENT SUPPLIED

| <i>Kit Components</i> | <i>State</i> | <i>Quantity</i> |
|--|--------------|-----------------|
| Antibody Coated Microtiter Strips | ready to use | 96 wells |
| Biotin Labelled Antibody Conc. (100x) | concentrated | 0.13 ml |
| Streptavidin-HRP Conjugate | ready to use | 13 ml |
| Master Standard | lyophilized | 2 vials |
| Quality Control HIGH | lyophilized | 2 vials |
| Quality Control LOW | lyophilized | 2 vials |
| Dilution Buffer | ready to use | 20 ml |
| Biotin-Ab Diluent | ready to use | 13 ml |
| Wash Solution Conc. (10x) | concentrated | 100 ml |
| Substrate Solution | ready to use | 13 ml |
| Stop Solution | ready to use | 13 ml |
| Product Data Sheet + Certificate of Analysis | - | 1 pc |

8. MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized (distilled) water
- Test tubes for diluting samples
- Glassware (graduated cylinder and bottle) for Wash Solution (Dilution Buffer)
- Precision pipettes to deliver 10-1000 μl with disposable tips
- Multichannel pipette to deliver 100 μl with disposable tips
- Absorbent material (e.g. paper towels) for blotting the microtiter plate after washing
- Vortex mixer
- Orbital microplate shaker capable of approximately 300 rpm
- Microplate washer (optional). [Manual washing is possible but not preferable.]
- Microplate reader with 450 ± 10 nm filter, preferably with reference wavelength 630 nm (alternatively another one from the interval 550-650 nm)
- Software package facilitating data generation and analysis (optional)

9. PREPARATION OF REAGENTS

- All reagents need to be brought to room temperature prior to use
- Always prepare only the appropriate quantity of reagents for your test
- Do not use components after the expiration date marked on their label

- Assay reagents supplied ready to use:

Antibody Coated Microtiter Strips

Stability and storage:

Return the unused strips to the provided aluminium zip-sealed bag with desiccant and seal carefully. Remaining Microtiter Strips are stable 3 months when stored at 2-8°C and protected from the moisture.

Streptavidin-HRP Conjugate

Dilution Buffer

Biotin-Ab Diluent

Substrate Solution

Stop Solution

Stability and storage:

Opened reagents are stable 3 months when stored at 2-8°C.

- **Assay reagents supplied concentrated or lyophilized:**

Human Vaspin Master Standard

Refer to the Certificate of Analysis for current volume of Dilution Buffer needed for reconstitution of standard!!!

Reconstitute the lyophilized Master Standard with Dilution Buffer just prior to the assay. Let it dissolve at least 15 minutes with occasional gentle shaking (not to foam). The resulting concentration of the human vaspin in the stock solution is **2 ng/ml**.

Prepare set of standards using Dilution Buffer as follows:

| <i>Volume of Standard</i> | <i>Dilution Buffer</i> | <i>Concentration</i> |
|---------------------------|------------------------|----------------------|
| Stock | - | 2 ng/ml |
| 250 µl of stock | 250 µl | 1 ng/ml |
| 250 µl of 1 ng/ml | 250 µl | 0.5 ng/ml |
| 250 µl of 0.5 ng/ml | 250 µl | 0.25 ng/ml |
| 250 µl of 0.25 ng/ml | 250 µl | 0.125 ng/ml |
| 250 µl of 0.125 ng/ml | 250 µl | 0.063 ng/ml |
| 250 µl of 0.063 ng/ml | 250 µl | 0.031 ng/ml |

Prepared Standards are ready to use, do not dilute them.

Stability and storage:

Do not store the diluted Standard solutions.

Quality Controls HIGH, LOW

Refer to the Certificate of Analysis for current volume of Dilution Buffer needed for reconstitution and for current Quality Control concentration!!!

Reconstitute each Quality Control (HIGH and LOW) with Dilution Buffer just prior to the assay. Let it dissolve at least 15 minutes with occasional gentle shaking (not to foam).

Reconstituted Quality Controls are ready to use, do not dilute them.

Stability and storage:

Do not store the reconstituted Quality Controls.

Biotin Labelled Antibody Conc. (100x)

Prepare the working Biotin Labelled Antibody solution by adding 1 part Biotin Labelled Antibody Concentrate (100x) with 99 parts Biotin-Ab Diluent. Example: 10 µl of Biotin Labelled Antibody Concentrate (100x) + 990 µl of Biotin-Ab Diluent for 1 strip (8 wells).

Stability and storage:

Opened Biotin Labelled Antibody Concentrate (100x) is stable 3 months when stored at 2-8°C.

Do not store the diluted Biotin Labelled Antibody solution.

Wash Solution Conc. (10x)

Dilute Wash Solution Concentrate (10x) ten-fold in distilled water to prepare a 1x working solution. Example: 100 ml of Wash Solution Concentrate (10x) + 900 ml of distilled water for use of all 96-wells.

Stability and storage:

The diluted Wash Solution is stable 1 month when stored at 2-8°C. Opened Wash Solution Concentrate (10x) is stable 3 months when stored at 2-8°C.

10. PREPARATION OF SAMPLES

The kit measures human vaspin in serum and heparin plasma.

Samples should be assayed immediately after collection or should be stored at -20°C. Mix thoroughly thawed samples just prior to the assay and avoid repeated freeze/thaw cycles, which may cause erroneous results. Avoid using hemolyzed or lipemic samples.

Dilute samples 3x with Dilution Buffer just prior to the assay, e.g. 50 µl of sample + 100 µl of Dilution Buffer for singlets, or preferably 100 µl of sample + 200 µl of Dilution Buffer for duplicates. **Mix well** (not to foam). Vortex is recommended.

Stability and storage:

Samples should be stored at -20°C, or preferably at -70°C for long-term storage. Avoid repeated freeze/thaw cycles.

Do not store the diluted samples.

See Chapter 13 for stability of serum samples when stored at 2-8°C, effect of freezing/thawing and effect of sample matrix (serum/plasma) on the concentration of human vaspin.

Note: It is recommended to use a precision pipette and a careful technique to perform the dilution in order to get precise results.

11. ASSAY PROCEDURE

1. Pipet **100 µl** of Standards, Quality Controls, Dilution Buffer (=Blank) and diluted samples, preferably in duplicates, into the appropriate wells. See *Figure 1* for example of work sheet.
2. Incubate the plate at room temperature (ca. 25°C) for **1 hour**, shaking at ca. 300 rpm on an orbital microplate shaker.
3. Wash the wells 5-times with Wash Solution (0.35 ml per well). After final wash, invert and tap the plate strongly against paper towel.
4. Add **100 µl** of Biotin Labelled Antibody solution into each well.
5. Incubate the plate at room temperature (ca. 25°C) for **1 hour**, shaking at ca. 300 rpm on an orbital microplate shaker.
6. Wash the wells 5-times with Wash Solution (0.35 ml per well). After final wash, invert and tap the plate strongly against paper towel.
7. Add **100 µl** of Streptavidin-HRP Conjugate into each well.
8. Incubate the plate at room temperature (ca. 25°C) for **30 min**, shaking at ca. 300 rpm on an orbital microplate shaker.
9. Wash the wells 5-times with Wash Solution (0.35 ml per well). After final wash, invert and tap the plate strongly against paper towel.
10. Add **100 µl** of Substrate Solution into each well. Avoid exposing the microtiter plate to direct sunlight. Covering the plate with e.g. aluminium foil is recommended.
11. Incubate the plate for **10 minutes** at room temperature. The incubation time may be extended [up to 20 minutes] if the reaction temperature is below than 20°C. Do not shake the plate during the incubation.
12. Stop the colour development by adding **100 µl** of Stop Solution.
13. Determine the absorbance of each well using a microplate reader set to 450 nm, preferably with the reference wavelength set to 630 nm (acceptable range: 550 - 650 nm). Subtract readings at 630 nm (550 - 650 nm) from the readings at 450 nm. **The absorbance should be read within 5 minutes following step 12.**

Note: If some samples and standard/s have absorbances above the upper limit of your microplate reader, perform a second reading at 405 nm. A new standard curve, constructed using the values measured at 405 nm, is used to determine vaspin concentration of off-scale standards and samples. The readings at 405 nm should not replace the readings for samples that were "in range" at 450 nm.

Note 2: Manual washing: Aspirate wells and pipet 0.35 ml Wash Solution into each well. Aspirate wells and repeat twice. After final wash, invert and tap the plate strongly against paper towel. Make certain that Wash Solution has been removed entirely.

| | strip 1+2 | strip 3+4 | strip 5+6 | strip 7+8 | strip 9+10 | strip 11+12 |
|----------|-----------------------|----------------|-----------|-----------|------------|-------------|
| A | Standard 2 | QC HIGH | Sample 7 | Sample 15 | Sample 23 | Sample 31 |
| B | Standard 1 | QC LOW | Sample 8 | Sample 16 | Sample 24 | Sample 32 |
| C | Standard 0.5 | Sample 1 | Sample 9 | Sample 17 | Sample 25 | Sample 33 |
| D | Standard 0.25 | Sample 2 | Sample 10 | Sample 18 | Sample 26 | Sample 34 |
| E | Standard 0.125 | Sample 3 | Sample 11 | Sample 19 | Sample 27 | Sample 35 |
| F | Standard 0.063 | Sample 4 | Sample 12 | Sample 20 | Sample 28 | Sample 36 |
| G | Standard 0.031 | Sample 5 | Sample 13 | Sample 21 | Sample 29 | Sample 37 |
| H | Blank | Sample 6 | Sample 14 | Sample 22 | Sample 30 | Sample 38 |

Figure 1: Example of a work sheet.

12. CALCULATIONS

Most microtiter plate readers perform automatic calculations of analyte concentration. The Standards curve is constructed by plotting the absorbance (Y) of Standards against the known concentration (X) of Standards, using the four-parameter algorithm. Results are reported as concentration of vaspin ng/ml in samples.

Alternatively, the *logit log* function can be used to linearize the standard curve, i.e. *logit* of the mean absorbance (Y) is plotted against *log* of the known concentration (X) of Standards.

The measured concentration of samples calculated from the standard curve must be multiplied by their respective dilution factor, because samples have been diluted prior to the assay, e.g. 0.15 ng/ml (from standard curve) x 3 (dilution factor) = 0.45 ng/ml.




Figure 2: Typical Standard Curve for Human Vaspin ELISA.

13. PERFORMANCE CHARACTERISTICS

➤➤ **Typical analytical data of BioVendor Human Vaspin ELISA are presented in this chapter**

- **Sensitivity**

Limit of detection (LOD) (defined as concentration of analyte giving absorbance higher than mean absorbance of blank* plus three standard deviations of the absorbance of blank: $A_{\text{blank}} + 3 \times \text{SD}_{\text{blank}}$) is calculated from the real human vaspin values in wells and is 0.01 ng/ml.

* Dilution Buffer is pipetted into Blank wells.

- **Limit of Assay**

Results exceeding human vaspin level of 2 ng/ml should be repeated with more diluted samples. Dilution factor needs to be taken into consideration in calculating the vaspin concentration.

- **Specificity**

The antibodies used in this ELISA are specific for human vaspin.

Sera of several mammalian species were measured in the assay. See results below.

For details please contact us at info@biovendor.com

| <i>Mammalian serum sample</i> | <i>Observed crossreactivity</i> |
|-------------------------------|---------------------------------|
| Bovine | no |
| Cat | no |
| Dog | no |
| Goat | no |
| Hamster | no |
| Horse | no |
| Monkey | no |
| Mouse | no |
| Pig | no |
| Rabbit | no |
| Rat | no |
| Sheep | no |

➤➤ **Presented results are multiplied by respective dilution factor**

• **Precision**

Intra-assay (Within-Run) (n=8)

| <i>Sample</i> | <i>Mean (ng/ml)</i> | <i>SD (ng/ml)</i> | <i>CV (%)</i> |
|---------------|-------------------------|-----------------------|-------------------|
| 1 | 0.613 | 0.039 | 6.5 |
| 2 | 0.210 | 0.018 | 8.7 |

Inter-assay (Run-to-Run) (n=6)

| <i>Sample</i> | <i>Mean (ng/ml)</i> | <i>SD (ng/ml)</i> | <i>CV (%)</i> |
|---------------|-------------------------|-----------------------|-------------------|
| 1 | 4.34 | 0.25 | 5.8 |
| 2 | 2.11 | 0.20 | 9.5 |

• **Spiking Recovery**

Serum samples were spiked with different amounts of human vaspin and assayed.

| <i>Sample</i> | <i>Observed (ng/ml)</i> | <i>Expected (ng/ml)</i> | <i>Recovery O/E (%)</i> |
|---------------|-----------------------------|-----------------------------|-----------------------------|
| 1 | 0.54 | - | - |
| | 0.76 | 0.92 | 82.6 |
| | 0.66 | 0.73 | 90.4 |
| | 0.62 | 0.63 | 98.4 |
| 2 | 0.29 | - | - |
| | 0.56 | 0.68 | 82.9 |
| | 0.45 | 0.49 | 91.8 |
| | 0.40 | 0.39 | 102.6 |

• **Linearity**

Serum samples were serially diluted with Dilution Buffer and assayed.

| <i>Sample</i> | <i>Dilution</i> | <i>Observed (ng/ml)</i> | <i>Expected (ng/ml)</i> | <i>Recovery O/E (%)</i> |
|---------------|-----------------|-----------------------------|-----------------------------|-----------------------------|
| 1 | - | 14.64 | - | - |
| | 2x | 7.08 | 7.32 | 96.7 |
| | 4x | 3.96 | 3.72 | 106.4 |
| | 8x | 2.04 | 1.80 | 113.3 |
| 2 | - | 1.47 | - | - |
| | 2x | 0.67 | 0.735 | 91.2 |
| | 4x | 0.35 | 0.370 | 94.6 |
| | 8x | 0.18 | 0.183 | 98.4 |

- **Effect of sample matrix**

EDTA, citrate and heparin plasmas were compared to respective serum samples from the same 10 individuals.

Results are shown below:

| Volunteer No. | Serum (ng/ml) | Plasma (ng/ml) | | |
|---|---------------|----------------|--------------|--------------|
| | | EDTA | Citrate | Heparin |
| 1 | 0.217 | 0.138 | 0.125 | 0.160 |
| 2 | 0.069 | 0.159 | 0.079 | 0.120 |
| 3 | 0.250 | 0.082 | 0.213 | 0.241 |
| 4 | 0.081 | 0.177 | 0.077 | 0.125 |
| 5 | 0.356 | 0.172 | 0.301 | 0.346 |
| 6 | 0.027 | 0.120 | 0.080 | 0.071 |
| 7 | 0.314 | 0.158 | 0.109 | 0.218 |
| 8 | 0.103 | 0.131 | 0.113 | 0.197 |
| 9 | 0.182 | 0.111 | 0.154 | 0.173 |
| 10 | 0.189 | 0.208 | 0.142 | 0.194 |
| Mean (ng/ml) | 0.179 | 0.145 | 0.139 | 0.185 |
| Mean Plasma/Serum (%) | - | 81.3 | 77.8 | 103.2 |
| Coefficient of Determination R² | - | 0.01 | 0.59 | 0.77 |




Figure 3: Vaspin levels measured using Human Vaspin ELISA from 10 individuals using serum, EDTA, citrate and heparin plasma, respectively.

- **Stability of samples stored at 2-8°C**

Samples should be stored at -20°C. However, no significant decline in concentration of human vaspin was observed in serum samples after 7 days when stored at 2-8°C. To avoid microbial contamination, samples were treated with ε-aminocaproic acid and sodium azide, resulting in the final concentration of 0.03% and 0.1%, respectively.

| <i>Sample</i> | <i>Incubation Temp, Period</i> | <i>Serum (ng/ml)</i> |
|---------------|--------------------------------|----------------------|
| 1 | -20°C | 0.429 |
| | 2-8°C, 1 day | 0.459 |
| | 2-8°C, 7 days | 0.456 |
| 2 | -20°C | 1.932 |
| | 2-8°C, 1 day | 1.683 |
| | 2-8°C, 7 days | 1.908 |
| 3 | -20°C | 0.261 |
| | 2-8°C, 1 day | 0.261 |
| | 2-8°C, 7 days | 0.168 |

- **Effect of Freezing/Thawing**

No significant decline was observed in concentration of human vaspin in serum samples after repeated (5x) freeze/thaw cycles. However it is recommended to avoid unnecessary repeated freezing/thawing of the samples.

| <i>Sample</i> | <i>Number of f/t cycles</i> | <i>Serum (ng/ml)</i> |
|---------------|-----------------------------|----------------------|
| 1 | 1x | 0.456 |
| | 3x | 0.435 |
| | 5x | 0.486 |
| 2 | 1x | 1.941 |
| | 3x | 1.926 |
| | 5x | 1.761 |
| 3 | 1x | 0.225 |
| | 3x | 0.258 |
| | 5x | 0.267 |

14. DEFINITION OF THE STANDARD

The recombinant human vaspin is used as the Standard. The human vaspin, produced in *E.coli*, is 46.7 kDa protein containing methionyl 408 amino acid residues of the human vaspin.

15. PRELIMINARY POPULATION AND CLINICAL DATA

The following results were obtained when serum samples from 128 unselected donors (82 men + 46 women) 15-82 years old were assayed with the Biovendor Human Vaspin ELISA in our laboratory:

- **Age and Sex dependent distribution of vaspin**

| Sex | Age (years) | n | Mean | SD | Min. | Max. | Median |
|-------|-------------|----|----------------|-------|-------|-------|--------|
| | | | Vaspin (ng/ml) | | | | |
| Men | 18 – 19 | 6 | 0.045 | 0.036 | 0 | 0.109 | 0.029 |
| | 20 – 49 | 60 | 1.698 | 5.713 | 0 | 38.15 | 0.071 |
| | 50 - 64 | 16 | 2.688 | 6.164 | 0.021 | 19.95 | 0.148 |
| Women | 15 – 18 | 2 | 0.459 | 0.42 | 0.039 | 0.879 | 0.459 |
| | 20 – 48 | 31 | 1.707 | 5.968 | 0.030 | 33.95 | 0.159 |
| | 50 - 82 | 13 | 0.139 | 0.071 | 0.039 | 0.336 | 0.117 |




Figure 4: Human vaspin concentration plotted against donor age and sex.

- **Reference range**

The data quoted in these instructions should be used for guidance only. It is recommended that each laboratory include its own panel of control sample in the assay. Each laboratory should establish its own normal and pathological reference ranges for vaspin levels with the assay.

16. METHOD COMPARISON

The Biovendor Human Vaspin ELISA was compared to another commercial ELISA immunoassay, by measuring 28 serum samples. The following correlation graph was obtained.




Figure 5: Method comparison.

17. TROUBLESHOOTING AND FAQs

»» Weak signal in all wells

Possible explanations:

- Omission of a reagent or a step
- Improper preparation or storage of a reagent
- Assay performed before reagents were allowed to come to room temperature
- Improper wavelength when reading absorbance

»» High signal and background in all wells

Possible explanations:

- Improper or inadequate washing
- Overdeveloping; incubation time with Substrate Solution should be decreased before addition of Stop Solution
- Incubation temperature over 30°C

»» High coefficient of variation (CV)

Possible explanation:

- Improper or inadequate washing
- Improper mixing Standards, Quality Controls or samples









18. REFERENCES

»» References to human vaspin:


- Li Q, Chen R, Moriya J, Yamakawa J, Sumino H, Kanda T and Takahashi T: A novel adipocytokine, visceral adipose tissue-derived serine protease inhibitor (Vaspin), and obesity. *J Int Med Res.* 2008; 36: 625-629
- Tan KB, Heutling D, Chen J, Farhatullah S, Adya R, Keay SD, Kennedy CR, Lehnert H and Randeve HS: Metformin decreases the adipokine vaspin in overweight women with polycystic ovary syndrome concomitant with improvement in insulin sensitivity and a decrease in insulin resistance. *Diabetes.* 2008; 57:1501-1507
- Youn B-S, Klötting N, Kratzsch J, Lee N, Park JW, Song E-S, Ruschke K, Oberbach A, Fasshauer M, Stumvoll M and Blüher M: Serum vaspin concentrations in human obesity and type 2 diabetes. *Diabetes.* 2008; 57:372-377
- Seeger J, Ziegelmeier M, Bachmann A, Lössner U, Kratzsch J, Blüher M, Stumvoll M and Fasshauer M: Serum levels of the adipokine vaspin in relation to metabolic and renal parameters. *J Clin Endocrinol Metab.* 2008; 93:247-251
- Wada J: Vaspin: a novel serpin with insulin sensitizing effects. *Expert Opin Investig Drugs.* 2008; 17(3): 327-333
- Hida K, Wada J, Eguchi J, Zhang H, Baba M, Seida A, Hashimoto I, Okada T, Yasuhara A, Nakatsuka A, Shikata K, Hourai S, Futami J, Watanabe E, Matsuki Y, Hiramatsu R, Akagi S, Makino H and Kanwar YS: Visceral adipose tissue-derived serine protease inhibitor: A unique insulin-sensitizing adipocytokine in obesity. *PNAS.* 2005; 102 (30): 10610-10615
- Klötting N, Berndt J, Kralisch S, Kovacs P, Fasshauer M, Schön MR, Stumvoll M and Blüher M: Vaspin gene expression in human adipose tissue: Association with obesity and type 2 diabetes. *BBRC.* 2005; 339: 430-436

»» For more references on this product see our WebPages at www.biovendor.com

19. EXPLANATION OF SYMBOLS

| | |
|---|---------------------------------------|
|  | Catalogue number |
|  | Content |
|  | Lot number |
|  | See instructions for use |
|  | Biological hazard |
|  | Expiry date |
|  | Storage conditions |
|  | Identification of packaging materials |


Assay Procedure Summary



| | | | | | | | | |
|----|----------|----------|----------|----------|----------|----------|----------|----------|
| | | | | | | | | |
| 12 | | | | | | | | |
| 11 | | | | | | | | |
| 10 | | | | | | | | |
| 9 | | | | | | | | |
| 8 | | | | | | | | |
| 7 | | | | | | | | |
| 6 | | | | | | | | |
| 5 | | | | | | | | |
| 4 | | | | | | | | |
| 3 | | | | | | | | |
| 2 | | | | | | | | |
| 1 | | | | | | | | |
| | A | B | C | D | E | F | G | H |

BioVendor

Research
and Diagnostic Products



HEADQUARTERS:
BioVendor - Laboratorní
medicína a.s.

Karasek 1767/1

621 00 Brno
CZECH REPUBLIC

Phone: +420-549-124-185
Fax: +420-549-211-460

E-mail: info@biovendor.com
sales@biovendor.com
Web: www.biovendor.com

EUROPEAN UNION:
BioVendor GmbH

Im Neuenheimer Feld 583

D-69120 Heidelberg
GERMANY

Phone: +49-6221-433-9100
Fax: +49-6221-433-9111

E-mail: infoEU@biovendor.com

USA, CANADA AND MEXICO:
BioVendor LLC

128 Bingham Rd.
Suite 1300

Asheville, NC 28806
USA

Phone: +1-828-575-9250
+1-800-404-7807
Fax: +1-828-575-9251

E-mail: infoUSA@biovendor.com

CHINA - Hong Kong Office:
BioVendor Laboratories Ltd

Room 4008
Hong Kong Plaza, No.188

Connaught Road West
Hong Kong, CHINA

Phone: +852-2803-0523
Fax: +852-2803-0525

E-mail: infoHK@biovendor.com

CHINA – Mainland Office:
BioVendor Laboratories Ltd

Room 2917, 29/F
R & F Ying Feng Plaza, No.2
Huaqiang road

Pearl River New Town
Guang Zhou, CHINA

Phone: +86-20-38065519
Fax: +86-20-38065529

E-mail: infoCN@biovendor.com