

CHEMICAL HAZARDS

Avoid contact with reagents containing TMB, hydrogen peroxide and sulfuric acid. If contacted with any of these reagents, wash with plenty of water. TMB is a suspected carcinogen.

SPECIMEN COLLECTION AND STORAGE

Serum: Approximately 0.2 mL of serum is required per duplicate determination. Collect 4-5 mL of blood into an appropriately labelled tube and allow it to clot. Centrifuge and carefully remove the serum layer. Store at 4°C for up to 24 hours or at -10°C or lower for longer time.

Plasma: Approximately 0.2 mL of plasma is required per duplicate determination. Collect 4-5 mL of blood into EDTA plasma tubes. Store at 4°C for up to 24 hours or at -10°C or lower for longer time.

Consider all human specimens as possible biohazardous materials and take precautions when handling.

SPECIMEN PREPARATION

Serum and plasma samples are loaded directly to the microplate wells; no specimen pretreatment is necessary.

REAGENTS AND EQUIPMENT NEEDED BUT NOT PROVIDED

1. Precision pipettes to dispense 25, 50, 100, 150 and 350 µL
2. Disposable pipette tips
3. Distilled or deionized water
4. Incubator set to 37°C.
5. Microplate well reader with a filter set at 450nm and an upper OD limit of 3.0 or greater

REAGENTS PROVIDED

1. Anti-Reverse T3 Polyclonal Antibody Coated Microplate -Break Apart Wells – Ready to Use	
Contents:	One 96 well (12x8) polyclonal antibody-coated microplate in a re-sealable pouch with desiccant.
Storage:	Refrigerate at 2-8°C
Stability:	12 months or as indicated on label

2. Reverse T3-Biotin Conjugate – Ready to Use	
Contents:	Reverse T3-Biotin conjugate in a protein-based buffer with a non-mercury preservative.
Volume:	13 mL/bottle
Storage:	Refrigerate at 2-8°C
Stability:	12 months in unopened bottle or as indicated on label.

3. Streptavidin-Horse Radish Peroxidase (HRP) Conjugate – Ready to Use	
Contents:	Streptavidin-HRP conjugate in a protein-based buffer with a non-mercury preservative.
Volume:	20 mL/bottle
Storage:	Refrigerate at 2-8°C
Stability:	12 months in unopened bottle or as indicated on label.

4. Reverse T3 Calibrators – Ready to Use	
Contents:	Six vials containing rT3 in a protein-based buffer with a non-mercury preservative. Prepared by spiking buffer with rT3 to the concentrations in labels. Typical calibrator concentrations*: 0; 0.02; 0.1; 0.4; 1 and 2 ng/mL. * Approximate value — please refer to vial labels for exact concentrations.
Volume:	Calibrators A-F: 1 mL/vial
Storage:	Refrigerate at 2-8°C
Stability:	12 months in unopened vials or as indicated on label


5. Reverse T3 Controls – Ready to Use	
Contents:	Two vials containing rT3 in a protein-based buffer with a non-mercury preservative. Prepared by spiking buffer with rT3 to the target concentration in QC certificate. Refer to vial labels for acceptable ranges.
Volume:	1 mL/vial
Storage:	Refrigerate at 2-8°C
Stability:	12 months in unopened vials or as indicated on label

6. Wash Buffer Concentrate – X10	
Contents:	One bottle containing buffer with a non-ionic detergent and a non-mercury preservative.
Volume:	50 mL/bottle
Storage:	Refrigerate at 2-8°C
Stability:	12 months or as indicated on label.
Preparation:	Dilute the wash buffer concentrate 1:10 in distilled or deionized water to prepare the <u>working wash buffer</u> . If one whole plate is to be used dilute 50 mL of the wash buffer concentrate in 450 mL of water.

7. TMB Substrate – Ready to Use	
Contents:	One bottle containing tetramethylbenzidine and hydrogen peroxide in buffer.
Volume:	16 mL/bottle
Storage:	Refrigerate at 2-8°C
Stability:	12 months or as indicated on label.

8. Stopping Solution – Ready to Use	
Contents:	One bottle containing 1M sulfuric acid.
Volume:	6 mL/bottle
Storage:	Refrigerate at 2-8°C
Stability:	12 months or as indicated on label.

SYMBOLS



ASSAY PROCEDURE

All reagents must reach room temperature before use. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.

1.	After all kit components have reached room temperature, mix gently by inversion. Prepare the working wash buffer (see wash buffer concentrate under the section "Reagents Provided").
2.	Remove the required number of strips from the microplate and assemble into a plate frame. Reseal the bag and return any unused strips to the refrigerator.
3.	Pipette 25 µL of each calibrator, control and specimen sample (serum or plasma) into correspondingly labelled wells in duplicate.
4.	Pipette 100 µL of the Reverse T3-Biotin conjugate into each well (the use of a multichannel pipette is recommended). Gently shake the microplate by hand for ten seconds to ensure complete mixing of the conjugate solution with the calibrators, controls and samples.
5.	Incubate the plate at 37°C for 1 hour. Do not cover the microplate.
6.	Wash the wells with 350 µL/well of working wash buffer solution 3 times . After washings tap the plate firmly against absorbent paper to remove any residual liquid (the use of an automatic strip washer is strongly recommended). —The accuracy of this assay depends on the correct execution of the washing procedure—
7.	Pipette 150 µL of the Streptavidin-HRP conjugate into each well. (the use of a multichannel pipette is recommended).
8.	Incubate the plate at 37°C for 30 minutes. Do not cover the microplate.
9.	Wash the wells 3 times using the same procedure as stated in step 6.
10.	Pipette 150 µL of the TMB substrate into each well at timed intervals (the use of a multichannel pipette is recommended).
11.	Incubate at 37°C for 15 minutes. Do not cover the microplate.
12.	Pipette 50 µL of stopping solution into each well at the same timed intervals as in step 10. (the use of a multichannel pipette is recommended). Gently shake the microplate by hand for ten seconds to ensure complete mixing of the stopping solution in the wells.
13.	Measure the absorbance at 450 nm with a microplate reader, within 20 minutes after addition of the stopping solution.

CALCULATIONS

- Calculate the mean optical density of each calibrator, control and specimen sample duplicate.
- Use a 4-parameter or 5-parameter curve with immunoassay software to generate the control and sample concentration results or draw a calibration curve on semi-log paper with the mean optical densities on the Y-axis and the calibrator concentrations on the X-axis and read the concentration of controls and samples off the calibrator curve.
- If a sample reads greater than 2 ng/mL report the result as ">2 ng/mL".
- To convert from ng/mL to ng/dL multiply the result by 100; to convert to nmol/L, multiply the ng/dL result by 0.01536 or the ng/mL result by 1.536.


TYPICAL TABULATED DATA

Sample data only. Do not use to calculate results

Calibrator	rT3 (ng/mL)	Mean OD (450 nm)
A	0	2.594
B	0.02	2.140
C	0.1	1.597
D	0.4	1.113
E	1	0.834
F	2	0.678
Unknown	0.15	1.458

TYPICAL CALIBRATION CURVE

Sample curve only. Do not use to calculate results.



PERFORMANCE CHARACTERISTICS

SENSITIVITY

The limit of detection (LoD) was determined from the analysis of 60 samples of the blank and a low value sample in two independent experiments and it was calculated as follows:

$LoD = \mu_B + 1.645\sigma_B + 1.645\sigma_S$, where σ_B and σ_S are the standard deviation of the blank and low value sample and μ_B is the mean value of the blank.

The Limit of Detection (LoD) was determined to be **0.009 ng/mL**.

SPECIFICITY

CROSS-REACTIVITY

The following compounds were tested for cross-reactivity with rT3 cross-reacting at 100%:

Steroid	%Cross Reactivity
rT3	100
T3	<0.001
T4	0.005
T2	0.004

INTERFERENT SUBSTANCES

The following substances did not show significant interference with the assay: hemoglobin up to 2 g/L, free and conjugated bilirubin up to 200 mg/L, triglycerides up to 5.5 mg/mL and Biotin up to 40 µg/mL.

INTRA-ASSAY PRECISION

Four serum samples were assayed 24 times each on the same calibrator curve. The results are tabulated below:

Sample	Mean (ng/mL)	SD (ng/mL)	%CV
1	0.089	0.0024	2.7
2	0.250	0.020	8.0
3	0.455	0.019	4.2
4	1.018	0.140	13.8

INTER-ASSAY PRECISION

Four serum samples were assayed in 20 different tests in the span of ten days. The results are tabulated below:

Sample	Mean (ng/mL)	SD (ng/mL)	CV%
1	0.127	0.016	12.6
2	0.304	0.038	12.5
3	0.469	0.057	12.2
4	0.847	0.083	9.8

RECOVERY

Spiked samples were prepared by adding defined amounts of rT3 to three patient serum samples. The results are tabulated below:

Sample	Observed Result (ng/mL)	Expected Result (ng/mL)	Recovery %
Serum Sample 1	0.437	-	-
+0.15 ng/mL	0.586	0.587	99.8
+0.30 ng/mL	0.728	0.737	98.8
+0.45 ng/mL	0.854	0.887	96.3
Serum Sample 2	0.327	-	-
+0.15 ng/mL	0.435	0.477	91.2
+0.30 ng/mL	0.596	0.627	95.1
+0.45 ng/mL	0.777	0.777	100.0
Serum Sample 3	0.149	-	-
+0.15 ng/mL	0.287	0.299	96.0
+0.30 ng/mL	0.460	0.449	102.4
+0.45 ng/mL	0.614	0.599	102.5

COMPARATIVE STUDIES

The DBC rT3 ELISA kit (y) was compared with the leading competing technology: Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS) (x). The comparison of 40 serum samples yielded the following linear regression results:

$$y = 0.79x + 0.01, \quad r = 0.96$$

REFERENCE VALUES

Reference values were obtained from commercial human specimens and calculated using a non-parametric method. Each laboratory must establish the range of reference values for their own population.

Group	n	Median (ng/mL)	95% Confidence Range (ng/mL)	Total range (ng/mL)
Serum adults	160	0.13	0.08-0.31	0.06-0.76
Plasma adults	120	0.14	0.08-0.29	0.049-0.65

LITERATURE


- Annewieke W et al. Thyroid Hormone Concentrations, Disease, Physical Function and Mortality in Elderly Men. *The Journal of Clinical Endocrinology and Metabolism*. **2005**, 90:6403-6409.
- Holtorf K *Thyroid Hormone Transport into Cellular Tissue*. *Journal of Restorative medicine*. **2014**, 3:53-68.
- Holtorf K *Peripheral Thyroid Hormone Conversion and Its Impact on TSH and Metabolic Activity*. *Journal of Restorative medicine*. **2014**, 3:30-52.
- Senese R et al. *Thyroid: Biological Actions of "Non-classical" Thyroid Hormones*, **2014**, 221:R1-R12
- Warner MH et al. *Mechanisms Behind the Non-thyroidal Illness Syndrome: An Update*. *Journal of Endocrinology*, **2010**, 205:1-13.
- Peeters RP et al. *Tissue Thyroid Hormone Levels in Critical Illness*. *Journal of Clinical Endocrinology and Metabolism*. **2005**, 90:6498-6507.
- Friberg L. et al. *Association Between Increased Levels of Reverse Triiodothyronine and Mortality after Acute Myocardial Infarction*. *American Journal of Medicine*, **2001**, 111:699-703.
- Pimentel et al. *Reverse T3 as a Parameter of Myocardial Function Impairment in Heart Failure*. *International Journal of Cardiology* **2010**, 145:52-53.
- Economidou F et al. *Thyroid Function During Critical Illness*. *Hormones* **2011**, 10:117-124.
- <https://www.nahypothyroidism.org/thyroid-hormone-transport/#reverseT3>

CONTACT INFORMATION

GENERAL:

 Diagnostics Biochem Canada (DBC) Inc.
41 Byron Ave,
Dorchester, Ontario, Canada, N0L 1G2
Tel (519) 268-8872
Fax (519) 268-7167
e-mail:dbc@dbc-labs.com
www.dbc-labs.com
An ISO 9001 and ISO 13485 Registered Company

IN EUROPE:

 Emergo Europe
Molenstraat 15
2513 BH, The Hague
The Netherlands